



# SKH040(1)

## LAUNDRY TROLLEY

Registration code: SKM-B



**SAIKANG** / Medical

Healthier  
your life

# SKH040(1)

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Size (LxWxH) 800×600×840 mm

### Technical configuration

3" Castor 4pcs



#### Structure

Φ25, Φ13mm stainless steel round tube, simple structure and practical; stainless steel material meets the needs of the vast market.



#### Material

One dirty clothes bag, waterproof cloth, waterproof and oil-proof, large capacity.



#### Dimension

Size: 800×600×840mm, Collects and transports dirty clothes, the body is light and easy to move.



#### Base | Steel

High-strength metal structure base frame, high bearing capacity and strong stability. Baosteel from Fortune Global 500.



#### Handle

Pushing handle, easy to move the trolley.



#### Castor

3" Bulging pole wheel, two with brakes, easy to install, stable and smooth.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue