



SKB039(C)

AMBULANCE

STRETCHER TROLLEY

Registration code: SKC-3B4



SAIKANG Medical

Healthier
your life

SKB039(C)

FEATURES AND HIGHLIGHTS




Technical parameters

▪ Full size (LxWxH)	1920x570x920mm
▪ Lowest size (LxWxH)	1920x570x280mm
▪ Back-rest lifting	0-60°
▪ Net weight	36kg
▪ Load bearing	≤135kg

Safety Lock

With safety lock, it must be locked up when transport the patients, and opened up when get on the ambulance.



Caster

5" silence caster. Advanced protection universal brake wheel, easy to move, simple to operate, reliable performance



Function

The mechanical foldaway legs are controlled by the handles on both ends of the stretcher, can be easily drop off when get on ambulance.



Fixing Devices

Can be locked by the fixing devices when on the ambulance car.



Structure

Premium alluminium alloy structure, more stable, but light weight



Leather / Safety Belts

Artificial leather, bright appearance, soft texture, abrasion resistance, folding resistance, acid and alkali resistance
With 2 safety belts, make patient safer



Back rest lifting

Back is adjustable which makes the wounded feel comfortable. Back rest lifting 0-60°





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue