



SKB3C01

ELECTRIC STAIR STRETCHER

Registration code: SKW-B2



SAIKANG/Medical

Healthier
your life

SKB3C01

FEATURES AND HIGHLIGHTS



Technical parameters

▪ Expanded size:	960x630x1150 mm
▪ Folding size:	880x630x485 mm
▪ Packing size:	1130x700x360 mm
▪ Capacity:	≤181kg
▪ Net weight:	42kg
▪ Gross weight:	52kg
▪ Motor:	200W
▪ Battery:	18AH
▪ Charger:	29.4V, 3A

Technical configuration

☒ Linear Motor	1 pc
☒ Lithium battery	1 pc
☒ Vest belt	1 set
☒ Wrench	1 pc
☒ Sponge mattress	2 pcs
☒ 8" wheel	2 pcs
☒ 12" wheel	2 pcs

Advantages

- ★ Powder coated steel material, the crawler is placed underneath, which makes climbing stairs more stable and labor-saving.
- ★ 24V, 200W brushless motor
- ★ Large-capacity lithium battery, long battery life.
- ★ Large-tooth shape crawler, stronger grip.
- ★ Honeycomb wheel, with shock absorption function, 8 inches front wheel is and 12 inches rear wheel.
- ★ Folding backrest.



Design

Scientific design, ergonomics, simple structure, simple operation, reliable and secure.

Widely used in the home for the old people, school security, hospitals, hotels etc for emergency evacuation.

**Material**

High strength powder coated steel construction provides durability and long life, foldable chair to storage and carry.

**Cushion**

PU sponge, detachable mattress. Height adjustable pillow. Foldable backrest.

**Armrest**

PU armrests, relieve hand pressure, easy to clean.

**Safety**

Circuit is inside the pipe, waterproof, noiseless during operation. Compliance with EU safety standards.

**Electric control system**

Built-in electronic speed control limits track speed during descent to a safe, manageable level. There is easy operation when climbing or descending stairs with intuitive, easy-to-use-style controls.

**Motor**

24V, 200W brushless motor, with cover, power system with soft start protection, make battery and motor using time longer, lower noise.

**Battery**

With detachable large capacity lithium battery, under seat, in line with international standards, not close to the human body, making users feel more comfortable.

**Function**

Glide tracked stair stretcher provides maximum safety when going up and down stairs using an electronically controlled motor and direct drive transmission with chair tracks.

**Safety belt**

Adjustable seat belts, to provide more protection for passengers.

**Front Wheels**

8" front wheel, rubber universal wheel, high quality and reliable, Honeycomb wheel, with shock absorption function.

**Rear Wheels**

12" rear wheel, rubber universal wheel, high quality and reliable, Honeycomb wheel, with shock absorption function.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 sales@chinasaikang.com
 www.saikangmedical.com



Company video



General catalogue