



# SKH005

## FAN-SHAPED INSTRUMENT TROLLEY

Registration code: SKM-B



**SAIKANG** Medical


Healthier  
your life

# SKH005

## FEATURES AND HIGHLIGHTS

### Technical parameters

|                    |                 |
|--------------------|-----------------|
| ▪ Size: L: (LxWxH) | 1550×460×900 mm |
| ▪ M: (LxWxH)       | 1400×460×900 mm |
| ▪ S: (LxWxH)       | 1200×460×900 mm |




### Technical configuration

☒ 3" Castor 5 pcs


#### Structural design

Φ25、Φ8mm stainless steel round tube, 1.2mm thick stainless steel board. The whole frame used high quality S.S. with the mirror surface polishing.



#### Application

Stainless steel instrument trolley, Multifunctional nursing trolley, Medical trolley.



#### Countertop structure

Two-layer countertop, high-quality stainless steel, equipped with three-sided stainless steel fence to prevent items from falling down.



#### Craftsmanship

Imported intelligent robot is used for welding forming, high precision, no blind spots, smooth surface, no burrs, beautiful and strong, high bearing capacity.



#### Castor

3 inch covered castors (diagonal brake) 5pcs, flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.





CONSTANT  
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE  
BRANCHES

#### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

#### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

#### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

#### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

#### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SAIKANG** Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ sales@chinasaikang.com  
🌐 www.saikangmedical.com



Company video



General catalogue