



SKH046-11

OVERBED TABLE




SAIKANG Medical

Healthier
your life

SKH046-11

FEATURES AND HIGHLIGHTS



Technical parameters

- External size (LxWxH) 800×395×855-1155mm

Technical configuration

- ☒ Air Spring for Height Adjustment 1pc
- ☒ Phenolic Resin Table Top 1pc

Material

Powder coated steel base, powder coated steel column.



Table Top

Table top is made of phenolic resin laminate with PVC cover.



Height Adjustment

Height adjustment by gas spring, used single-handed according to the convenience of the user.



Casters

Easy to install, more stable and firmer
High quality medical caster, 1.5 inch /2 inch mute wheel.



Matching Design

It can be used with the bedside table, which is easy to disassemble. The combination of dining table and cabinet functions can serve two purposes and save space.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue