



SKR-J10

PATIENT MONITOR CART

Registration code: SKM-A



SAIKANG Medical

Healthier
your life

SKR-J10

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 505×430×950 mm

Technical configuration

☒ Big Basket	1pc
☒ Tray	1pc
☒ Castor	4pcs



Column

Aluminum alloy column, the top equipment mesa.



Computer Holder

Concise and compact for 13" patient monitor, Proposal Table load ≤ 15 kg.



Big Basket

With basket on the back, easy to place goods, Environmentally plastic is non-toxic and easy to clean.



Tray

Built-in ABS tray, ABS tray size: 388mm (L) × 238mm (W) × 50mm (H), the height of the ABSI tray above the ground is 579.8mm.



Base

Base is ABS covered high quality metal structure, featured with both stability and easy-maintenance.



Castor

4pcs of 75mm double side casters, each with brake.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue