



# SKR-MT625

## MEDICINE TROLLEY

Registration code: SKM-B



**SAIKANG** Medical

Healthier  
your life

# SKR-MT625

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Size (LxWxH)	625×475×930 mm
----------------	----------------

### Technical configuration

<input checked="" type="checkbox"/> Dust Baskets	1pc
<input checked="" type="checkbox"/> Utility Container	1pc
<input checked="" type="checkbox"/> Central Lock	1set
<input checked="" type="checkbox"/> Trash Can	2pcs
<input checked="" type="checkbox"/> Castor	4pcs



#### Drawers

1-4 layer small drawer, fifth floor middle drawer, sixth floor large drawer, can put different sizes of medicine boxes and other medical materials, and equipped with central control lock.



#### Trash Can

There are two ABS dirt buckets of the car body, different colors can distinguish different kinds of dirt.



#### Heavy Duty Construction

Wheels assembled on the same vertical line as trolley four columns, more stable.



#### Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



#### Powder Coated Column

Make the structure of the trolley more stable and firm.



#### Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley.



#### File Box

Side attached file box, for storage of paper work.



#### Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.



#### Drawer Card Holder

card insert from downside, easy keep clean.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SAIKANG** Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎ +86 18021231901 · 🕒 24/7  
✉ sales@chinasaikang.com  
🌐 www.saikangmedical.com



Company video



General catalogue