



R000w

ONE FUNCTION MANUAL BED

Registration code: SK-CD1



SAIKANG Medical

Healthier
your life

R000w

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxWxH)	2000x1000x490 mm
▪ Mattress platform	1950x950 mm
▪ Safe working load	200 kg
▪ Back-rest adjustment	0-75° (±10°)



Technical configuration

☒ Powder coated Bed Ends	1set
☒ IV Pole	1pc
☒ IV Pole Prevision	1p

Optional

☐ 5" Noiseless Castors	
☐ Integrated detachable guardrail	



Steels | Painting


Baosteel from Fortune Global500.
11 Process epoxy painting, ASTM
testing anti-bacterial, paint thickness
0.12 mm, brightness 60°, paint can
resist 50kg impact.

**Welding**

Panasonic robotic ensure, 360° full
smooth welding

**Bed Surface**

Reticulated bed platform.

**Head & Foot Board**

Powder coated head&foot board,
simple installation,Removable.

**Handle**


Equipped with a handle to manually
control the angle of back lift.

**IV Pole**

Full stainless steel IV pole, bearing
15 kg.

**Optional**

Integrated detachable guardrail.
Model: G7A

**Optional**

5" noiseless castors, two with break.
Model: C4B





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue