



X01-4

BABY CRIB

Registration code: SK-AD3

SAIKANG Medical

Healthier
your life

X01-4

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Bed Height Min.	840 mm
▪ Bed Height Max.	960 mm
▪ Reversed Trendelenburg	0-10°
▪ Safe working load	25 kg

Technical configuration

<input checked="" type="checkbox"/> Hydraulic Pump	1pc
<input checked="" type="checkbox"/> Bassinet	1pc
<input checked="" type="checkbox"/> 4" Castors	4pcs

Design

Modern design, ergonomics and provide a comfortable environment for babies.



Baby Basin

The transparent baby basin is detachable, and the environmentally friendly plastic is non-toxic and odorless, also can be disinfected and sterilized at a temperature above 100°C.



Wheel

4-inch casters, wear-resistant, non-slip, can rotate 360° freely.



Reversed trendelenburg

It is convenient for communication between mother and baby, and also provides a suitable angle for baby feeding.




Weighing system

The baby's weight can be measured at any time, and the indoor temperature can be measured at the same time.



Lifting system

Adjust the height through the hydraulic system, stable and smooth.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue