



A99-8

GYNECOLOGICAL EXAM COUCH

Registration code: SKL-B






SAIKANG Medical

Healthier
your life

A99-8

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External Size (LxWxH)	1450x600x580-1000 mm
▪ Back Plate Up/Down	70°/20°
▪ Seat board upward	40°
▪ Main Voltage AC	220V±22V,50Hz±1Hz
▪ Patient Weight Capacity	185kg

Technical configuration

☒ Pull Handle	2pcs
☒ Knee Curtch	2pcs
☒ Pedal Remote Controller	1pc
☒ Stainless Steel Basin	1pc
☒ Instrution Manual	1pc

Lowest Position of Table Top

580mm



Highest Position of Table Top

1000mm



Back Plate Upward

0°-70°



Back Plate Downward

0°-20°



Seat Board Upward

0°-40°



Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.

**Function**


Double using leg support , can be used as boatstretcher.

**Pedal Remote Controller**

Famous brand motor,Various actions are controlled by remote controller.

**Pull Handle**

The hand bar can be dismantled and movable.

**Stainless Steel Basin**


Equipped with stainless steel sewage basin and other various accessories.

**Mattress**

The mattress adopts high density memory sponge and is once formed,with no seam.

**Wheel**

Each wheel equipped with single brake which is easy for moving obstetric table.

**Knee Curtch**



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue