



# SKH032

## TREATMENT TROLLEY

Registration code: SKM-B



**SAIKANG** Medical

Healthier  
your life

# SKH032

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 670×450×1700 mm

### Technical configuration

- ☒ 3" Castor 4pcs



#### Application

Stainless steel infusion treatment trolley; Multifunctional nursing trolley; Medical Trolley.



#### Structural design

The structure is compact and practical, use high-quality stainless steel, safe and non-toxic, easy to clean, smooth surface, not easy to corrode, and long service life.



#### Countertop structure

Two-layer countertop and drawer structure is adopted. The countertop is high-quality stainless steel, equipped with three-sided stainless steel fence to prevent items from falling down.



#### Basin

The product adopts a dirt basin design, the dirt can be put into the basin, it is convenient to disassemble and assemble, and the basin can be taken out directly.



#### Drawer

The drawer is equipped with high-quality silent slide rails, which can be drawn smoothly, effectively solves the problem of easily sliding out or blocking during the drawing process.



#### Castors

3 inch covered castors (diagonal brake), flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎ +86 18021231901 · 🕒 24/7  
✉ sales@chinasaikang.com  
🌐 www.saikangmedical.com



Company video



General catalogue