

Good night,
better days

SK-D07-4

ELECTRIC HOMECARE BED

User's safety and comfort from long term care bed

www.saikangmedical.com



SAIKANG Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
www.saikangmedical.com



POSTURE FUNCTION



USER'S SAFETY AND COMFORT FROM LONG TERM CARE BED

The physical health and psychological well-being are hot topics for elderly people. Safety is an important quality feature for a nursing bed. Beds' systematic fall prevention can reduce the risk of falling out of the bed.

Our innovative solutions can reduce falls as well as the consequences associated with falls. It is a wonderful match of safety, mobility, height adjustable and coziness. All these features cater to the needs of people who are spend a lot time in bed. The adjustability makes beds more appropriate than a standard bed for providing comfortable support and enabling caregivers to deliver quality care.



Safe working load

200_{kg}



Back-rest adjustment

65°



Leg lifting adjustment

25°



Knee-rest adjustment

45°



Right-turning adjustment

50°



Left-turning adjustment

50°

ELECTRIC + MANUAL

Bed functions can be driven by cranks or motors. Then, the bed can work in the case of power failure.



Motor and controller

Adopted Jiechang motors, brings transmission technology, meet 3C certificates. You can set up and use by the remote controller surprisingly simple.



Rocker

Extra rocker for emergency use when the power cut out.



Guardrail and end of bed

Paired with imported PP side rails, product quality stable and reliable. Buffer lifting, anti pinch, provides the patient with ultimate safety when using this bed.



Safe working load


200_{kg}

Bedboard Design

Its durable structure comes with a weight capacity of 200kg, this bed set accommodates most patients. An open architecture frame makes cleaning the bed quick and easy.

Bed Material

Double bed frame construction, welded with high quality steel, steel hinge connection.





Caster

Medical silent casters with brake system,ensure smooth mobility of the bed in any direction.



Colors

This home care Bed features a wooden color guardrails and head foot board which will compliment any home setting.



Pooholes

Hidden poohole design for one-handed push and pull quickly. The size of the hole is 200*270mm, designed by ergonomics.



POSTURE FUNCTION


Good night, better days

Built with incredible attention to detail, this home care bed is versatile enough to comfortably accommodate any aspect of day-to-day life.



Back&Knee-rest Lifting

This function is perfect for anyone aided by the bed helping to stand and reducing the pressure of knee.



Back-rest Lifting

Back-rest Lifting will bring you to a seated position for for watching TV, eating a meal, or visiting with family and provide a moderate height for caregivers.



Knee-rest Lifting

This unique function gives the patient support and promoting blood circulation and help with leg pain and swelling.

RESET

One-button reset

Press the one-button reset button after each power failure to prevent the bed damage while turning.



ACCESSORIES


Mattress

Removable cotton cover is durable, Built-in 4cm brown mat + 2cm latex, water resistant, anti-microbial.



Automatic turn over function

Through pressing the button to control turning over Interval of 30 minutes or 45 minutes. Press the button again can stop this function.



Left-turning&right-turning

Helps users turning easily, improving blood circulation and reducing the chance of skin shearing or pressure sores.



Dining table (optional)

Moving table easily for comfort while reading, writing or eating.



Infusion rack

The infusion rack is made of stainless steel, and its height can be adjustable.



Hair washing basin

Hair washing basin with drain pipe, draining the dirty water easily for users after used.



Bedpan

The bedpan is made of raw material by injection molding, so it has the advantage of convenient cleaning.

Good night, better days

SAIKANG Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
www.saikangmedical.com

+86 18021231901

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

export@saikangmedical.com

www.saikangmedical.com



Product video (scan / press QR code)

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.