



# A3008A-1

## MANUAL OPERATING TABLE

Registration code: SKF-C




**SAIKANG** Medical

Healthier  
your life

# A3008A-1

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ External Size (LxWxH):	2100x480x750-950 mm
▪ Built-in Kidney Bridge:	0-100 mm
▪ Back Plate Up/Down:	75°/5°
▪ Head Plate Up/Down:	30°/90°
▪ Trendelenburg/Reverse Trendelenburg:	30°/30°
▪ Lateral Tilt Left/Right:	20°/20°
▪ Leg Plate Down:	90°
▪ Patient Weight Capacity:	250 kg
▪ X-rays translucent:	optional with extra charge
▪ C-arm compatible:	NO
▪ Multi-functional orthopedic frame:	YES

### Technical configuration

<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Long Fixing Clamp	1pair
<input checked="" type="checkbox"/> Fixing Clamp	8pcs

Lowest Position of Table Top  
**770mm**



Highest Position of Table Top  
**1020mm**




Back Plate Upward  
**0°-80°**



Back Plate Downward  
**0°-10°**



Trendelenburg  
Reversed Trendelenburg  
**25°/25°**



Lateral Tilt Left  
**0°-20°**



Lateral Tilt Right  
**0°-20°**



Leg Plate Down  
**0°-90°**



**Application**

This operating table is a comprehensive multi-function type, suitable for medical units to perform various operations, all movements are controlled by the manual cranks.

**Material**

The base, column and side hoder are all made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.

**Control System**


Foot treadle brake device which is easily for moving operating table.

**Built-In Kidney Bridge**





Built-In kidney bridge, more convenient and nice appearance.

**Hydraulic**

The height adjustment controlled by hydraulic pump.

**Mattress**

PU cover, custom mattress, easy disassembly and cleaning, thickness of mattress: 6cm

**Shoulder Support with Cushion****Anesthesia Screen Frame****Arm Support with Cushion****Leg Support with Cushion****Body Support with Cushion**



## **CONSTANT INNOVATION**



Design



Research



Manufacturing



Maintenance



Service

## **WORLDWIDE BRANCHES**

### South America Branch

#### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

#### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

#### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

#### Uk Healthward International Limited

DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

#### LLC 'SAIKANG'

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya  
Naberezhnaya, Russia  
TEL: +7 (919)300-00-06

### Mid East Branch


#### Saikang Medical Supplies (Mid East) LLC

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

 +86 18021231901 ·  24/7

 sales@saikangmedical.com

 www.saikangmedical.com



Company video



General catalogue